Sodio ialuronato alto peso e niacinamide per uso intra-articolare

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Italcert srl
Viale Sarca, 336 - 20126 MILANO (MI)
Country: Italy

Hyaluronic acid sodium salt HW and niacinamide for intra-articular use

PRODUCT DESCRIPTION

Playal is a sterile, biodegradable, isotonic, injectable gel, for intra-articular use. *Playal* consists of High chain (1.8 -2.6 x 10⁶ Dalton) hyaluronic acid, obtained from *Streptococcus equi* bacteria, formulated to a concentration of 20 mg/ml in a physiologic buffer. *Playal* is characterized by viscoelastic properties, therefore allows to facilitate the normalization of the viscosity of the synovial fluid present in the intra-articular cavity.

Each box contains one syringe of *Playal* and a product leaflet. A set of two labels showing the batch number is contained in the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

In addition a implant card and an implant card label (if needed) is given with the package. Niacinamide is an excipient used to protect hyaluronic acid and to give it more stability.

COMPOSITION

2 mL syringes: Sodium hyaluronate 20 mg /ml, sodium chloride, niacinamide, sodium dihydrogen phosphate dihydrate, dibasic sodium phosphate dodecahydrate, WFI grade water.

EXPECTED CLINICAL BENEFITS

Primary endpoint of clinical benefits is based on the percentage of patients with therapeutic success, defined as alleviation and/or remission of pain symptoms and difficulty to perform daily normal actions for those who are prone to joint degenerative diseases or trauma.

PERFORMANCE REQUIREMENTS

Playal is based on high molecular weight hyaluronic acid and niacinamide which is an excipient for protection of hyaluronic acid. Hyaluronic acid in his form (high viscosity solution with viscoelastic and lubricant properties) has the capability to adsorb the impact and to avoid additional shock in already traumatized joints. Niacinamide, as an excipient, has the purpose to protect hyaluronic acid and to guarantee its durability to achieve its purpose. The protection operate by niacinamide is in terms of oxidation and in terms of protection from the action of hyaluronidase. In this way hyaluronic acid could achieve more durability. (In vitro test: E.IT.076 2019/1536 and 2019/1537).

INTENDED PURPOSE

Playal is a synovial fluid substitute which, thanks to its viscoelastic and lubricant properties, promotes the restoration of rheological conditions of the joints, altered in degenerative or post-traumatic conditions. The product, improving the characteristics of the synovial fluid, exerts a protective action of the joints and helps the improvement of joint function and the reduction of pain symptoms. *Playal* acts only at the joint where it is injected, without exerting any systemic action.

POPULATION OF USE

Population of use: 18-85, male and female. Normal people and professional in sport.

WARNINGS - PRECAUTIONS FOR USE

Playal is suitable only for intra-articular injections and must only be dispensed by a doctor who has received specific training on the intra-articular injections technique.

The structures in which the products could be administered are the following:

- orthopedics clinics;
- Hospital
- rehabilitation clinics in which there is one (or more) specialized doctor/doctors in intra-articular injections;
- private clinics;
- health center;
- rheumatologists' clinics in which there is one (or more) specialized doctor/doctors in intra-articular injections.

Before use, check the integrity of the syringe and the expiration date. Do not use needles other than those listed.

The product should not be injected in the presence of a infected or severely inflamed joint. The infiltration must be avoided in the case of infections in place or inflammatory conditions of the skin in proximity of the injection.

As no clinical experience is available for the use of Hyaluronan in children, treatment with *Playal* is not recommended in these cases. Better not use during concomitant therapies or other therapies with other medical device and other medicinal product in light of the fact that no evidence is been collected regarding concomitant therapies with the medical device or with product which could be related to the metabolism of this medical device. The device must not be used in pregnant or breastfeeding women, as has not been tested in such cases. therapies. The device must not be used in patients with ascertained individual hypersensitivity to any of the ingredients of the product.

Playal must be administered with caution to patients with a history of hypersensitivity to other drugs and patients with liver failure or impairment.

After the intra-articular injection it is advisable to recommend to the patient to avoid physical activities demanding stress for the articulation and resume normal activities after a few days.

After the intra-articular injection it is advisable to recommend to the patient to avoid driving for at least 48 hours.

Playal is a disposable product, the quality and sterility are guaranteed only if the syringe is sealed. Any residue must be discarded and not reused even after new sterilization. Do not use the product if the package is already opened or damaged. The assembled syringe must be discarded immediately after use, regardless of whether or not the solution has been completely administered. After use, dispose according to applicable national practice.

DO NOT RE-USE: quality and sterility can be guarantee only for an originally closed syringes. The re- use of the product creates a potential infection risk for patients and users.

DO NOT RE-STERILIZED

If the sterile package will be unintentionally opened before the use, discharge the syringes.

There are no known risks of overdose because the product is well tolerated but it is better not to exceed a total amount of 16 mL of gel injection during the year to avoid accumulation in joints.

INCOMPATIBILITIES

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as solutions of benzalkonium chloride. Contact between *Playal* and these substances should be therefore avoided.

SIDE EFFECTS

There may be some temporary side reactions following injection of *Playal* such as pain, stiffness, warmth, redness or swelling. These secondary manifestations may be relieved by applying ice on the treated articulation. Usually these effects disappear after a short time (3- 4 days maximum). If symptoms persist, seek medical attention. Any other unwanted side effects associated with the *Playal* injection must be reported to the doctor.

As for any intra-articular treatment, septic arthritis may rarely occur when general precautions for injections are not observed or the site of injection is not aseptic. More rarely, cases of hematoma and arthralgia are reported, as well as skin rash, such as urticaria and pruritus.


METHODS OF USE

Remove any joint effusion before injecting *Playal*; for the removal of the effusion and the injection of *Playal* the same needle must be used. Remove the protective cap of the syringe, with particular attention to avoid contact with the opening. Firmly screw the needle, of diameter between 18 and 22 g, at the collar of Luer lock, following the instruction given below. Before injection the site should be treated with appropriate disinfectant. Inject *Playal* adopting aseptic technique. Inject only into the joint cavity. In case of necessity use imaging technique to make the injections. It is recommended to perform an initial cycle of 3 weekly treatments in case of severe condition. Possibly followed by maintenance sessions, according to the medical prescription. Based on the condition of the patients the cycle could be different: 2 injections every 2 weeks followed, if necessary by maintenance session after 6 months in case of less severe condropathy.


The dosage of the syringes could be chosen based on the joint's cavity space and the severity of condition. For syringe of 4 mL (one shot) one injections is recommended followed by maintenance session if necessary.

INSTRUCTIONS FOR ASSEMBLY OF THE SYRINGE NEEDLE

A. Carefully unscrew the cap of the tip of the syringe, being particularly careful to avoid contact with the opening.

B. Gently grip the needle guard and mount the needle on the luer-lok® mount, screwing it tight until a slight counter-pressure is felt in order to ensure an airtight grip and prevent leakage of the liquid during administration. The cover of the needle must be held firmly during the procedure to avoid any leaking and malfunctioning. The cover must be pull and not unscrew to avoid the loss of connection between the needle and the syringes.

If the luer is not well connected or there is lack of product of the aspect is not in compliance the syringes must be discarded.

When the injections is about to be prepared the operator must not take the syringes from the plunger rod. This is to avoid a disassembling.

When the operator is about to make the injection he cannot touch the needle or the final part of the luer to avoid to contaminate the area.

After the use, the syringe and the needle must be throw away immediately without touching the needle to avoid injuries.

ADVERSE EVENTS

Unknown, if the correct methods of use and instructions given in the illustrative are observed. Any use outside the normal use described in this leaflet should be avoided.

Any adverse event or serious accident must be reported to the Manufacturer in the contacts below and to the Competent Authority of the member state in which the patient / user is established.

LANGUAGES AVAILABLE ON THE LEAFLET

The leaflet contains different languages.

The leaflet will be available in addition on the website: www.innate.it

STORAGE

Store *Playal* at 2–25 °C (36–77 °F) in a dry place in the original box. Protect from light, heat and frost. Keep out of reach of children.

CONTENTS OF THE PACK

The product is available in 2 different model:

- Pre-filled syringe containing 2 ml of non pyrogenic gel, sterilised using moist heat (*Playal* 40 with 0,8% mg of niacinamide) and implant card, implant label and leaflet
- Pre-filled syringe containing 4 ml of non pyrogenic gel, sterilised using moist heat (*Playal* 80 with 0,5% of niacinamide), implant card, implant label and leaflet

Summary of Safety and clinical performance could be found in Eudamed database or in the website of the manufacturer, in addition to a document for the important information for the patient.

HOW TO USE IMPLANT CARD

Implant card is given with the medical device.

The aim of introducing an IC has been to achieve three main objectives:

1. Enable the patient to identify the implanted devices and to get access to other information related to the implanted device (e.g. via EUDAMED, and other websites).
2. Enable patients to identify themselves as persons requiring special care in relevant situations e.g. security checks.
3. Enabling e.g. emergency clinical staff or first responder to be informed about special care/needs for relevant patients in case of emergency situations.

Some information is been already printed by the Manufacturer. The blank space must be filled by the final user in handwritten way as the following:

- Patient Name or patient ID.
The final user must be filled this space with the name of the patients (full name).
- Date of implantation.
The final user must put the date of the impantation of the medical device (in this case injection date).
- Name and Address of the implanting healthcare institution/provider.
Who perform the injections and in which hospital/clinic.

The other part of implant card could be already printed by the manufacturer or could be blank with the sentence "area for applying the IC label".

In the second case, the IC label is given with the package. This label must be applied by the doctor or the patient after the visit in the dedicated blank part.

The final version if implant card with the applied label must be kept by the patient.

LAST REVISED

Revision 04 of 08/11/2024

INDEX OF SYMBOLS

- Read instruction for use
- Limit of temperature
- Sterile fluid path sterilize agent steam/moist heat
- Don't use if package is damage
- Do no re-sterilize
- Keep dry
- Single use
- Manufacturer
- Manufacturing date
- Keep away from heating sources
- Batch number
- Expiration date
- Symbol for medical device
- Catalogue number
- Single sterile barrel system (sterilized with moist heat) with protective packaging outside
- CE mark
- Do not leave the packaging free in the environment
- Symbol for identification of UDI
- Website information of the Manufacturer

INNATE s.r.l. Viale Industria, 11-13 15067 Novi Ligure (AL) – ITALIA
info@innate.it – +39 0143 2645 – Lun-Ven 8:00-17:00

Italcert srl
Viale Sarca, 336 – 20126 MILANO (MI)
Country: Italy

CE 0426